

Anti-Ddx43 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nxz2

Item Number: CM0020724

Introduction

High-quality rabbit polyclonal antibody against DDX43 (CT13/HAGE), validated for Western blot and ELISA. Suitable for human, mouse, and rat samples. Detects a 73kDa ATP-dependent RNA helicase with cancer/testis antigen expression. Ideal for oncology research.

[Learn More](#)

Parameter	Value
Product Name	DDX43 Rabbit pAb
Remarks/Alias	CT13; HAGE; DDX43
Species	Human
GeneID (Human)	55510
GeneID	55510
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 429-648 of human DDX43 (NP_061135.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NXZ2
Protein Weight	73kDa
Shipping	Ice bag