


Anti-Ddx46 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q7L014

Item Number: CM0025114



Introduction

Rabbit polyclonal antibody against human DDX46 (PRP5/PRPF5), a spliceosomal RNA helicase. Validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for spliceosome assembly and pre-mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	DDX46 Rabbit pAb
Remarks/Alias	Prp5; PRPF5; DDX46
Species	Human
GeneID (Human)	9879
GeneID	9879
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L014
Protein Weight	117kDa
Shipping	Ice bag