


Anti-Ddx5 (P68) Polyclonal Antibody For Wb, Ihc-P, Elisa - P17844

Item Number: CM0014109



Introduction

Rabbit polyclonal antibody against human DDX5 (p68), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for RNA helicase and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	DDX5 Rabbit pAb
Remarks/Alias	p68, HLR1, G17P1, HUMP68, DDX5
Species	Human
GeneID (Human)	1655
GeneID	1655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-614 of human DDX5 (NP_004387.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17844
Protein Weight	69kDa
Shipping	Ice bag