


Anti-Ddx5 Polyclonal Antibody For Wb, Elisa - P17844

Item Number: CM0029430



Introduction

High-quality rabbit polyclonal antibody against human DDX5 (p68 RNA helicase), validated for WB and ELISA, with cross-reactivity to mouse. Ideal for studying transcriptional regulation, splicing, and circadian rhythms. SWISS: P17844.

[Learn More](#)

Parameter	Specification
Product Name	DDX5 Rabbit pAb
Remarks/Alias	p68; HLR1; G17P1; HUMP68; DDX5
Species	Human
Gene ID (Human)	1655
Gene ID	1655
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-614 of human DDX5 (NP_004387.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17844
Protein Weight	69kDa
Shipping	Ice bag