

Anti-Ddx5 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P17844

Item Number: CM0002099



Introduction

High-specificity rabbit monoclonal antibody targeting DDX5 (p68 helicase), validated for WB, IHC-P, IF/ICC, IP, ELISA with human, mouse, rat cross-reactivity. Ideal for IVD assay development, transcriptional regulation, and splicing research.

[Learn More](#)

Parameter	Value
Product Name	DDX5 Rabbit mAb
Remarks/Alias	p68; HLR1; G17P1; HUMP68; DDX5
Species	Human
Gene ID (Human)	1655
Gene ID	1655
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17844
Protein Weight	69kDa
Shipping	Ice bag