

Anti-Ddx56 Polyclonal Antibody For Wb, Elisa - Q9Ny93

Item Number: CM0025664

Introduction

Rabbit polyclonal antibody against human DDX56, validated for WB and ELISA. Suitable for studies on ribosome biogenesis, innate immunity, and viral infection research. Reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	DDX56 Rabbit pAb
Remarks/Alias	DDX21; DDX26; NOH61; DDX56
Species	Human
GeneID (Human)	54606
GeneID	54606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 318-547 of human DDX56 (NP_061955.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NY93
Protein Weight	62kDa
Shipping	Ice bag