


Anti-Ddx6 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26196

Item Number: CM0024127



Introduction

Rabbit polyclonal antibody against DDX6, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for studying P-bodies, mRNA decapping, and autophagy regulation.

[Learn More](#)

Parameter	Value
Product Name	DDX6 Rabbit pAb
Remarks/Alias	P54; RCK; HLR2; IDDILF; Rck/p54; DDX6
Species	Human
Gene ID (Human)	1656
Gene ID	1656
Immunogen	Synthetic peptide corresponding to amino acids 384-483aa of human DDX6 (NP_001244120.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26196
Protein Weight	54kDa
Shipping	Ice bag