

Anti-Decorin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P07585

Item Number: CM0006648

Introduction

Discover our specific rabbit monoclonal antibody against human Decorin, suitable for WB, IF-P, IHC-P, ELISA applications. Reacts with Human and Rat. Essential tool for studying extracellular matrix biology and proteoglycan function.

[Learn More](#)

Parameter	Specification
Product Name	Decorin Rabbit mAb
Remarks/Alias	CSCD; PG40; PGII; PGS2; DSPG2; SLRR1B
Species	Human
GeneID (Human)	1634
GeneID	1634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-359 of human Decorin (NP_001911.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P07585
Protein Weight	40kDa
Shipping	Ice bag