


Anti-Defa5/Defensin Alpha 5 Polyclonal Antibody For Wb, Elisa - Q01523

Item Number: CM0018950



Introduction

Rabbit polyclonal antibody targeting human DEFA5 (Defensin alpha 5), suitable for Western blotting and ELISA. Validated for human species, this antibody detects a ~10kDa antimicrobial peptide involved in innate immune defense.

[Learn More](#)

Parameter	Value
Product Name	DEFA5 Rabbit pAb
Remarks/Alias	DEF5; HD-5; DEFA5
Species	Human
GeneID (Human)	1670
GeneID	1670
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-94 of human DEFA5 (NP_066290.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q01523
Protein Weight	10kDa
Shipping	Ice bag