

Anti-Defb121 Rabbit Polyclonal Antibody For If-P, Ihc-P, Elisa - Q5J5C9

Item Number: CM0015744

Introduction

DEFB121 Rabbit Polyclonal Antibody validated for IF-P, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Targets beta-defensin 121, an antimicrobial peptide of the male reproductive tract. Ideal for reproductive immunology research.

[Learn More](#)

Parameter	Value
Product Name	DEFB121 Rabbit pAb
Remarks/Alias	DEFB21; ESC42RELC; DEFB121
Species	Human
Gene ID (Human)	245934
Gene ID	245934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human DEFB121 (NP_001011878.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5J5C9
Protein Weight	8kDa
Shipping	Ice bag