

Anti-Defb132 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q7Z7B7

Item Number: CM0012855

Introduction

Rabbit polyclonal antibody targeting human DEFB132 (beta-defensin 132), validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Recognizes a secreted antimicrobial peptide involved in innate immunity. Ideal for defensin research.

[Learn More](#)

Parameter	Value
Product Name	DEFB132 Rabbit pAb
Remarks/Alias	BD-32; DEFB32; HEL-75; UNQ827; DEFB-32; KFLL827; DEFB132
Species	Human
Gene ID	400830
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human DEFB132 (NP_997352.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z7B7
Protein Weight	11kDa
Shipping	Ice bag