

Anti-Defb132 Rabbit Polyclonal Antibody For Wb/Elisa - Q7Z7B7

Item Number: CM0017278

Introduction

Rabbit polyclonal antibody targeting human beta-defensin 132 (DEFB132, BD-32). Suitable for WB and ELISA; cross-reacts with mouse. 11 kDa protein. Ideal for antimicrobial peptide research.

[Learn More](#)

Parameter	Value
Product Name	DEFB132 Rabbit pAb
Remarks/Alias	BD-32; DEFB32; HEL-75; UNQ827; DEFB-32; KFLL827; DEFB132
Species	Human
Gene ID	400830
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q7Z7B7
Protein Weight	11kDa
Shipping	Ice bag