

Anti-Defb4A (Beta-Defensin 2) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - O15263

Item Number: CM0026452



Introduction

DEFB4A Rabbit Polyclonal Antibody for WB, IF-P, ELISA targets human beta-defensin 2 (BD-2), an antimicrobial peptide. Ideal for innate immunity and infection research.

[Learn More](#)

Parameter	Value
Product Name	DEFB4A Rabbit pAb
Remarks/Alias	BD-2, SAP1, DEFB2, DEFB4, HBD-2, DEFB-2, DEFB102
Species	Human
Gene ID (human)	1673
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-64 of human DEFB4A (NP_004933.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	O15263
Protein Weight	7kDa
Shipping	Ice bag