

Anti-Degs1 Rabbit Polyclonal Antibody For Wb, Elisa - 015121

Item Number: CM0013284

Introduction

This anti-DEGS1 rabbit polyclonal antibody detects human DEGS1, a sphingolipid delta(4)-desaturase crucial for converting sphinganine to sphingosine. Validated for Western blot and ELISA applications. Ideal for metabolic and neurological research.

[Learn More](#)

Parameter	Value
Product Name	DEGS1 Rabbit pAb
Remarks/Alias	MLD; DEGS; DES1; Des-1; FADS7; HLD18; MIG15; DEGS-1; DEGS1
Species	Human
GeneID (human)	8560
GeneID	8560
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-323 of human DEGS1 (NP_003667.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15121
Protein Weight	38kDa
Shipping	Ice bag