

Anti-Delta-Sarcoglycan (Sgcd) Rabbit Monoclonal Antibody For Wb, Elisa - Q92629

Item Number: CM0010258

Introduction

Recombinant rabbit monoclonal antibody against Delta-sarcoglycan (SGCD) validated for WB and ELISA. Detects human and mouse SGCD (~32 kDa). Ideal for muscular dystrophy and sarcolemma integrity studies.

[Learn More](#)

Parameter	Value
Product Name	Delta-sarcoglycan (SGCD) Rabbit mAb
Remarks/Alias	SGD; DAGD; 35DAG; CMD1L; SGCDP; LGMDR6; SG-delta; Delta-sarcoglycan (SGCD)
Species	Human
Gene ID (Human)	6444
Gene ID	6444
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 190-289 of human Delta-sarcoglycan (SGCD) (Q92629).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92629
Protein Weight	32kDa
Shipping	Ice bag