

Anti-Delta-Sarcoglycan (Sgcd) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92629

Item Number: CM0026927

Introduction

Rabbit polyclonal antibody targeting human delta-sarcoglycan (SGCD), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for sarcoglycan complex and muscular dystrophy studies.

[Learn More](#)

Parameter	Value
Product Name	Delta-sarcoglycan (SGCD) Rabbit pAb
Remarks/Alias	SGD; DAGD; 35DAG; CMD1L; SGCDP; LGMDR6; SG-delta; Delta-sarcoglycan (SGCD)
Species	Human
GeneID	6444
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-290 of human Delta-sarcoglycan (SGCD) (NP_000328.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92629
Protein Weight	32kDa
Shipping	Ice bag