

Anti-Denn (Madd) Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8Wxg6

Item Number: CM0002281

Introduction

Rabbit monoclonal antibody against human DENN (MADD), validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for synaptic trafficking, MAPK pathway, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	DENN Rabbit mAb
Remarks/Alias	DENN; IG20; DEEAH; RabGEF; NEDDISH; RAB3GEP
Species	Human
GeneID (human)	8567
GeneID	8567
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1548-1647 of human DENN (Q8WVG6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WVG6
Protein Weight	183kDa
Shipping	Ice bag