

Anti-Dennd4A Polyclonal Antibody For Wb And Elisa - Q7Z401

Item Number: CM0015754

Introduction

Rabbit polyclonal antibody against DENND4A (MYCPBP) for WB and ELISA. Cross-reacts with human, mouse. Targets ~209 kDa protein. Suitable for MYC promoter regulation and RAB10 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	DENND4A Rabbit pAb
Remarks/Alias	IRLB; MYCPBP; DENND4A
Species	Human
Gene ID (Human)	10260
Gene ID	10260
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1500-1650 of human DENND4A (NP_005839.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z401
Protein Weight	209kDa
Shipping	Ice bag