


Anti-Denr Polyclonal Antibody For Wb, Ihc-P, Elisa - O43583

Item Number: CM0015062



Introduction

High-quality rabbit polyclonal antibody against human DENR (DRP1/SMAP-3), validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for translation reinitiation studies.

[Learn More](#)

Parameter	Value
Product Name	DENR Rabbit pAb
Remarks/Alias	DRP; DRP1; SMAP-3; DENR
Species	Human
GeneID	8562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DENR (NP_003668.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43583
Protein Weight	22kDa
Shipping	Ice bag