


Anti-Dermcidin (Dcd) Rabbit Polyclonal Antibody For Wb, Elisa - P81605

Item Number: CM0026526



Introduction

High-quality rabbit polyclonal antibody targeting human Dermcidin (DCD), validated for WB and ELISA. Cross-reacts with mouse. Supports antimicrobial peptide and cancer biomarker research. Immunogen: aa 20-110. Protein weight: 11kDa.

[Learn More](#)

Parameter	Value
Product Name	DCD Rabbit pAb
Remarks/Alias	PIF; AIDD; DSEP; HCAP; DCD-1; DCD
Species	Human
GeneID (Human)	117159
GeneID	117159
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human DCD (NP_444513.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P81605
Protein Weight	11kDa
Shipping	Ice bag