

Anti-Desmin Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P17661

Item Number: CM0014100

Introduction

High-quality Desmin rabbit polyclonal antibody suitable for WB, IF-P, IHC-P, ELISA applications. Cross-reacts with human, mouse, and rat samples. Research-use IVD raw material for studying sarcomere integrity and intermediate filaments.

[Learn More](#)

Parameter	Value
Product Name	Desmin Rabbit pAb
Remarks/Alias	CSM1; CSM2; CDCD3; LGMD1D; LGMD1E; LGMD2R; Desmin
Species	Human
GeneID (human)	1674
GeneID	1674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human Desmin (NP_001918.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17661
Protein Weight	54kDa
Shipping	Ice bag