


Anti-Desmocollin 1 Monoclonal Antibody For Wb, Elisa - Q08554

Item Number: CM0007242



Introduction

Recombinant rabbit monoclonal antibody against Desmocollin 1, validated for WB and ELISA. Reacts with human and rat. Ideal for cell adhesion and desmosome research. SWISS: Q08554.

[Learn More](#)

Parameter	Value
Product Name	Desmocollin 1 Rabbit mAb
Remarks/Alias	CDHF1; DG2/DG3; Desmocollin 1
Species	Human
GeneID (human)	1823
GeneID	1823
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q08554
Protein Weight	100kDa
Shipping	Ice bag