

Anti-Desmocollin 3 Rabbit Monoclonal Antibody For Wb, Elisa - Q14574

Item Number: CM0003041

Introduction

Anti-Desmocollin 3 rabbit monoclonal antibody, raised against a synthetic peptide (aa 797-896 of human Q14574), validated for WB and ELISA. Suitable for studying desmosomal adhesion and epidermal barrier function.

[Learn More](#)

Parameter	Value
Product Name	Desmocollin 3 Rabbit mAb
Remarks/Alias	DSC; DSC1; DSC2; DSC4; CDHF3; HT-CP; Desmocollin 3
Species	Human
Gene ID (Human)	1825
Gene ID	1825
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 797-896 of human Desmocollin 3 (Q14574).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14574
Protein Weight	100kDa
Shipping	Ice bag