


Anti-Desmoglein 2 Rabbit Monoclonal Antibody For Wb, Elisa - Q14126

Item Number: CM0005152



Introduction

High-quality Anti-Desmoglein 2 Rabbit Monoclonal Antibody (Q14126) validated for WB and ELISA. Detects human Desmoglein 2, a desmosomal cadherin involved in cell adhesion and signaling. Ideal for cell junction and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Desmoglein 2 Rabbit mAb
Remarks/Alias	HDGC; CDHF5; Desmoglein 2
Species	Human
GeneID (Human)	1829
GeneID	1829
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1001-1118 of human Desmoglein 2 (Q14126)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14126
Protein Weight	122kDa
Shipping	Ice bag