

Anti-Desmoglein 3 Rabbit Polyclonal Antibody For Wb And Elisa - P32926

Item Number: CM0028789



Introduction

Rabbit polyclonal antibody against human Desmoglein 3 (DSG3, PVA). Validated for Western blot (WB) and ELISA. Suitable for studying desmosome biology, cell adhesion, and pemphigus vulgaris. Recombinant immunogen (aa 749-999).

[Learn More](#)

Parameter	Value
Product Name	Desmoglein 3 Rabbit pAb
Remarks/Alias	PVA; ABOLM; CDHF6; Desmoglein 3
Species	Human
GeneID (Human)	1830
GeneID	1830
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 749-999 of human Desmoglein 3 (NP_001935.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P32926
Protein Weight	108kDa
Shipping	Ice bag