

## Anti-Desmoglein 3/Pva Monoclonal Antibody For Wb, Ihc-P, Elisa - P32926

Item Number: CM0009213



### Introduction

Rabbit monoclonal antibody targeting human Desmoglein-3 (DSG3, PVA) amino acids 50-615, validated for WB, IHC-P, and ELISA. Research desmosome adhesion, pemphigus vulgaris autoantigen, and hair follicle biology.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Desmoglein 3/PVA Rabbit mAb                                                                                                  |
| Remarks/Alias    | PVA; ABOLM; CDHF6; DSG3                                                                                                      |
| Species          | Human                                                                                                                        |
| GeneID (human)   | 1830                                                                                                                         |
| GeneID           | 1830                                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 50-615 of human Desmoglein 3/PVA (NP_001935.2) |
| Source           | Rabbit                                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                                        |
| Application      | WB, IHC-P, ELISA                                                                                                             |
| Cross Reactivity | Human                                                                                                                        |
| SWISS            | P32926                                                                                                                       |
| Protein Weight   | 108kDa                                                                                                                       |
| Shipping         | Ice bag                                                                                                                      |