


Anti-Desmoplakin Polyclonal Antibody For Wb, Ihc-P, Elisa - P15924

Item Number: CM0025907



Introduction

Validate Desmoplakin expression with this rabbit polyclonal antibody. Suitable for WB, IHC-P, ELISA in human, mouse, rat. Research desmosome junctions, cell adhesion, and skin/heart biology. Target: P15924.

[Learn More](#)

Parameter	Value
Product Name	Desmoplakin Rabbit pAb
Remarks/Alias	DP; DCWHKTA; Desmoplakin
Species	Human
GeneID (human)	1832
GeneID	1832
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15924
Protein Weight	332kDa
Shipping	Ice bag