


Anti-Desmoplakin Polyclonal Antibody For Wb, Ihc-P, Elisa - P15924

Item Number: CM0015647



Introduction

Rabbit polyclonal antibody against Desmoplakin, validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat. Target weight 332kDa. Ideal for desmosome and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Desmoplakin Rabbit pAb
Remarks/Alias	DP; DCWHKTA; Desmoplakin
Species	Human
GeneID (Human)	1832
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2609-2830 of human Desmoplakin (NP_004406.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15924
Protein Weight	332kDa
Shipping	Ice bag