

## Anti-Det1 Polyclonal Antibody For Wb, Elisa - Q7L5Y6

Item Number: CM0028985



### Introduction

Rabbit polyclonal anti-DET1 antibody for WB and ELISA; cross-reactive with human, mouse, rat. DET1 is a component of the E3 ligase complex involved in JUN degradation. Ideal for ubiquitination and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | DET1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | DET1                                                                                                             |
| Species          | Human                                                                                                            |
| GeneID (human)   | 55070                                                                                                            |
| GeneID           | 55070                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 15-240 of human DET1 (NP_060466.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q7L5Y6                                                                                                           |
| Protein Weight   | 64kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |