

Anti-Dexi Polyclonal Antibody For Wb, Elisa - O95424

Item Number: CM0018264

Introduction

High-quality DEXI Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human DEXI (MYLE, 10 kDa) protein. Ideal for cardiac, pulmonary disease research and dexamethasone-induced expression studies.

[Learn More](#)

Parameter	Value
Product Name	DEXI Rabbit pAb
Remarks/Alias	MYLE; DEXI
Species	Human
GeneID (human)	28955
GeneID	28955
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human DEXI (NP_054734.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95424
Protein Weight	10kDa
Shipping	Ice bag