

Anti-Dffa Polyclonal Antibody For Wb, Elisa - 000273

Item Number: CM0016837



Introduction

DFFA (DFF-45/ICAD) rabbit polyclonal antibody for western blot and ELISA. Detects human and mouse DFFA, a key apoptosis inhibitor of caspase-activated DNase. Ideal for cell death research.

[Learn More](#)

Parameter	Value
Product Name	DFFA Rabbit pAb
Remarks/Alias	DFF1; ICAD; DFF-45; DFFA
Species	Human
Gene ID	1676
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-331 of human DFFA (NP_004392.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00273
Protein Weight	37kDa
Shipping	Ice bag