


Anti-Dffb Polyclonal Antibody For Wb, Ihc-P, Elisa - 076075

Item Number: CM0016788



Introduction

Rabbit polyclonal antibody against human DFFB (CAD/CPAN/DFF40), suitable for WB, IHC-P, and ELISA. Validated in human and mouse samples. Ideal for apoptosis and DNA fragmentation research.

[Learn More](#)

Parameter	Value
Product Name	DFFB Rabbit pAb
Remarks/Alias	CAD, CPAN, DFF2, DFF40, DFF-40, DFFB
Species	Human
GeneID (human)	1677
GeneID	1677
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 109-338 of human DFFB (NP_004393.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O76075
Protein Weight	39kDa
Shipping	Ice bag