

Anti-Dgat1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75907

Item Number: CM0030973



Introduction

Rabbit polyclonal antibody against human DGAT1, validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for studies of triacylglycerol synthesis and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	DGAT1 Rabbit pAb
Remarks/Alias	ARAT; DGAT; ARGPI; DIAR7; DGAT1
Species	Human
GeneID	8694
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 201-300 of human DGAT1 (NP_036211.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75907
Protein Weight	55kDa
Shipping	Ice bag