

Anti-Dgat2 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Pd7

Item Number: CM0015877

Introduction

This Anti-DGAT2 Rabbit Polyclonal Antibody is validated for Western blot and ELISA, reacting with human, mouse, and rat DGAT2. DGAT2 is the essential acyltransferase catalyzing the terminal step in triacylglycerol synthesis, making this antibody a valuable tool for lipid metabolism and obesity studies.

[Learn More](#)

Parameter	Value
Product Name	DGAT2 Rabbit pAb
Remarks/Alias	ARAT; HMFN1045; GS1999FULL; DGAT2
Species	Human
GeneID (Human)	84649
GeneID	84649
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 289-388 of human DGAT2 (NP_115953.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PD7
Protein Weight	44 kDa
Shipping	Ice bag