

Anti-Dgcr6L Polyclonal Antibody For Wb, Elisa - Q9By27

Item Number: CM0017981

Introduction

Rabbit polyclonal antibody against DGCR6L, suitable for Western blot and ELISA applications. Reacts with human, mouse, and rat samples. DGCR6L is a nuclear protein potentially involved in neural crest cell migration.

[Learn More](#)

Parameter	Value
Product Name	DGCR6L Rabbit pAb
Remarks/Alias	DGCR6; DGCR6L
Species	Human
GeneID (Human)	85359
GeneID	85359
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human DGCR6L (NP_150282.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY27
Protein Weight	25kDa
Shipping	Ice bag