


Anti-Dgcr8 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Wyq5

Item Number: CM0005325



Introduction

DGCR8 Rabbit mAb (Q8WYQ5) is a monoclonal antibody validated for WB and ELISA, recognizing human and mouse DGCR8. Ideal for miRNA biogenesis and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	DGCR8 Rabbit mAb
Remarks/Alias	Gy1; pasha; DGCRK6; C22orf12; DGCR8
Species	Human
GeneID (Human)	54487
GeneID	54487
Immunogen	Recombinant fusion protein containing amino acids 436-764 of human DGCR8 (Q8WYQ5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WYQ5
Protein Weight	86kDa
Shipping	Ice bag