

Anti-Dgka Rabbit Monoclonal Antibody For Wb, Elisa - P23743

Item Number: CM0003924



Introduction

Rabbit monoclonal DGKA antibody for WB and ELISA. Cross-reacts with human, mouse, rat. Targets diacylglycerol kinase alpha, a key enzyme in lipid signaling and DAG/PA balance.

[Learn More](#)

Parameter	Value
Product Name	DGKA Rabbit mAb
Remarks/Alias	DAGK; DAGK1; DGK-alpha; DGKA
Species	Human
Gene ID (Human)	1606
Gene ID	1606
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DGKA (NP_001336.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23743
Protein Weight	83kDa
Shipping	Ice bag