

Anti-Dgka Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P23743

Item Number: CM0003049



Introduction

Rabbit monoclonal antibody against human diacylglycerol kinase alpha (DGKA). Validated for Western blot, IHC paraffin, and ELISA applications. Useful for investigating lipid signaling, DAG/PA metabolism, and lymphocyte function.

[Learn More](#)

Parameter	Value
Product Name	DGKA Rabbit mAb
Remarks/Alias	DAGK; DAGK1; DGK-alpha; DGKA
Species	Human
GeneID (Human)	1606
GeneID	1606
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DGKA (NP_001336.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P23743
Protein Weight	83kDa
Shipping	Ice bag