

Anti-Dgka Rabbit Polyclonal Antibody For Wb, Elisa - P23743

Item Number: CM0028211



Introduction

DGKA rabbit polyclonal antibody for Western blot and ELISA, validated in human, mouse, and rat. Targets diacylglycerol kinase alpha, a key regulator of DAG/PA signaling and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	DGKA Rabbit pAb
Remarks/Alias	DAGK; DAGK1; DGK-alpha; DGKA
Species	Human
Gene ID (human)	1606
Gene ID	1606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DGKA (NP_958852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23743
Protein Weight	83kDa
Shipping	Ice bag