


Anti-Dgkz (Diacylglycerol Kinase Zeta) Rabbit Polyclonal Antibody For Wb And Elisa - Q13574

Item Number: CM0023429



Introduction

Rabbit polyclonal antibody targeting human diacylglycerol kinase zeta (DGKZ). Validated for WB and ELISA, with cross-reactivity to mouse. Ideal for lipid signaling, T cell activation, and actin regulation research.

[Learn More](#)

Parameter	Value
Product Name	DGKZ Rabbit pAb
Alias	DAGK5; DAGK6; DGK-ZETA; hDGKzeta; DGKZ
Species	Human
GeneID	8525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 808-1117 of human DGKZ (NP_001099010.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13574
Protein Weight	104kDa
Shipping	Ice bag