

Anti-Dglucy Monoclonal Antibody For Wb, Ihc-P, Elisa - Q7Z3D6

Item Number: CM0009400

Introduction

Recombinant rabbit monoclonal antibody against human DGLUCY (C14orf159), validated for WB, IHC-P, ELISA. Detects endogenous D-glutamate cyclase in human and mouse samples. Ideal for mitochondrial metabolism research.

[Learn More](#)

Parameter	Value
Product Name	DGLUCY Rabbit mAb
Remarks/Alias	C14orf159; DGLUCY
Species	Human
Gene ID (Human)	80017
Gene ID	80017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 269-616 of human DGLUCY (NP_001095836.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z3D6
Protein Weight	66kDa
Shipping	Ice bag