


Anti-Dguok Polyclonal Antibody For Wb, Elisa - Q16854

Item Number: CM0026076



Introduction

Rabbit polyclonal antibody targeting mitochondrial deoxyguanosine kinase (DGUOK) for WB and ELISA applications. Validated in human, mouse, and rat samples. Ideal for mitochondrial DNA synthesis and nucleoside metabolism research.

[Learn More](#)

Parameter	Value
Product Name	DGUOK Rabbit pAb
Remarks/Alias	dGK; NCPH; NCPH1; PEOB4; MTDPS3; DGUOK
Species	Human
Gene ID (Human)	1716
Gene ID	1716
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 40-277 of human DGUOK (NP_550438.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16854
Protein Weight	32kDa
Shipping	Ice bag