


Anti-Dhfr Polyclonal Antibody For Wb - P00374

Item Number: CM0018820



Introduction

DHFR rabbit polyclonal antibody detects human and mouse DHFR in WB, IF/ICC, and ELISA applications. Key for folate metabolism and thymidylate biosynthesis research.

[Learn More](#)

Specification	Details
Product Name	DHFR Rabbit pAb
Remarks/Alias	DYR, DHFR1, DHFRP1, DHFR
Species	Human
Gene ID (Human)	1719
Immunogen	Recombinant fusion protein containing amino acids 1-187 of human DHFR (NP_000782.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00374
Protein Weight	21kDa
Shipping	Ice bag