

Anti-Dhodh Polyclonal Antibody For Wb, Ihc-P, Elisa - Q02127

Item Number: CM0027835



Introduction

Rabbit polyclonal DHODH antibody validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat dihydroorotate dehydrogenase, a key mitochondrial enzyme in de novo pyrimidine biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	DHODH Rabbit pAb
Remarks/Alias	URA1; POADS; DHodehase; DHODH
Species	Human
GeneID (human)	1723
GeneID	1723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 116-395 of human DHODH (NP_001352.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02127
Protein Weight	43kDa
Shipping	Ice bag