

Anti-Dhodh Rabbit Polyclonal Antibody For Wb, Elisa - 035435

Item Number: CM0029020

Introduction

Rabbit polyclonal antibody targeting dihydroorotate dehydrogenase (DHODH) for WB and ELISA applications. Cross-reacts with human and mouse, supporting studies in pyrimidine metabolism and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	DHODH Rabbit pAb
Remarks/Alias	2810417D19Rik
Species	Mouse
GeneID (human)	1723
GeneID	56749
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-395 of mouse DHODH (NP_064430.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O35435
Protein Weight	43 kDa
Shipping	Ice bag