


Anti-Dhps Rabbit Polyclonal Antibody For Wb And Elisa - P49366

Item Number: CM0029787



Introduction

DHPS Rabbit Polyclonal Antibody validated for WB and ELISA against human, mouse, and rat DHPS (41kDa). Deoxyhypusine synthase catalyzes the essential hypusine modification of eIF-5A for translation regulation.

[Learn More](#)

Parameter	Value
Product Name	DHPS Rabbit pAb
Remarks/Alias	DS; DHS; MIG13; NEDSSWI; DHPS
Species	Human
GeneID (Human)	1725
GeneID	1725
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-369 of human DHPS (NP_001921.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49366
Protein Weight	41kDa
Shipping	Ice bag