

Anti-Dhrs1 Polyclonal Antibody For Wb And Elisa - Q96Lj7

Item Number: CM0022213



Introduction

Rabbit polyclonal antibody targeting human DHRS1 (SDR19C1) for Western blot and ELISA. Cross-reacts with mouse. Ideal for investigating steroid and xenobiotic metabolism and endoplasmic reticulum function.

[Learn More](#)

Parameter	Value
Product Name	DHRS1 Rabbit pAb
Remarks/Alias	SDR19C1; DHRS1
Species	Human
Gene ID	115817
Immunogen	Recombinant protein corresponding to amino acids 1-262 of human DHRS1 (NP_001129522.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96LJ7
Protein Weight	34kDa
Shipping	Ice bag