

## Anti-Dhrs11 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Uwp2

Item Number: CM0015142



### Introduction

Rabbit polyclonal antibody targeting human DHRS11 (SDR24C1). Validated for Western blot and ELISA, with cross-reactivity to mouse. Useful for steroid hormone metabolism and androgen biosynthesis research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | DHRS11 Rabbit pAb                                                                                                  |
| Remarks/Alias    | ARPG836; SDR24C1; spDHRS11; DHRS11                                                                                 |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 79154                                                                                                              |
| GeneID           | 79154                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DHRS11 (NP_077284.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | Q6UWP2                                                                                                             |
| Protein Weight   | 28kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |