

Anti-Dhrs2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13268

Item Number: CM0028581

Introduction

Rabbit polyclonal antibody targeting human DHRS2 (HEP27/SDR25C1). Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Investigates dicarbonyl reduction, p53 stabilization, and cancer biology.

[Learn More](#)

Parameter	Specification
Product Name	DHRS2 Rabbit pAb
Remarks/Alias	HEP27; SDR25C1; DHRS2
Species	Human
GeneID (Human)	10202
GeneID	10202
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13268
Protein Weight	30kDa
Shipping	Ice bag