

Anti-Dhrs9 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bpw9

Item Number: CM0030577

Introduction

DHRS9 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Detects a 35kDa protein. Ideal for retinol and steroid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	DHRS9 Rabbit pAb
Remarks/Alias	RDHL; RDH15; RDHTBE; SDR9C4; RDH-TBE; RETSDR8; hRoDH-E2; 3ALPHA-HSD; 3-alpha-HSD; DHRS9
Species	Human
GeneID (human)	10170
GeneID	10170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-319 of human DHRS9 (NP_001135743.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BPW9
Protein Weight	35kDa
Shipping	Ice bag